



*LeadCross*TM Guide Wire

LeoMed
乐奥医疗 为生命健康服务


GI Endoscopic Accessories Series 

LeadCross™ Guide Wire


Used for guiding the catheter or other devices to the digestive tract for interventional diagnostic and surgery endoscopically.




- **Soft tip**
Easy to insert and less tissue damage.



- **Good Visibility**
Colorful coating is easy to identify. Convenient for following the device movement under endoscope.



- Unique super-elastic nitinol core designed to promote excellent kink resistance, Dynamic nitinol core intended to promote pushability for easy advancement.



- Hydrophilic coating designed to facilitate crossing of tight lesions.

Product Specification

REF	Diameter (inch)	Working length (mm)	Tip shape	Hardness
L023520SS00	0.035	2000	Straight	Super stiff
L023526SS00	0.035	2600		
L023545SS00	0.035	4500		
L023520SN00	0.035	2000	Straight	Normal stiff
L023526SN00	0.035	2600		
L023545SN00	0.035	4500		