



# Single-use Balloon Dilatation Catheter Wire-guided Progressive Dilatation Balloon



- **Excellent Deliverability**
  - Soft and tapered tip allows less invasive advancement across the stricture.
  - Low-profile design reduces resistance and increases crossability and trackability.



- **Operational Efficiency**
  - Facilitate placement within tight strictures and tortuous anatomy.
  - Optimized balloon deflation time.



- **Optimized Balloon Design**
  - The balloon component is designed to offer stability during stricture dilation.
  - Semi-compliant material engineered to adapt to biliary duct.

## Specifications

| REF      | REF         |            |                        | REF         |            |             | Minimum Working Channel (mm) | Unit/Box |
|----------|-------------|------------|------------------------|-------------|------------|-------------|------------------------------|----------|
|          | Length (mm) | Diam. (mm) | Nominal Pressure (atm) | Length (mm) | Diam. (mm) | Structure   |                              |          |
| KZ0655G2 | 55          | 6-7-8      | 3-4.5-6                | 1800        | 2.4        | Wire-guided | 2.8                          | 1        |
| KZ0855G2 | 55          | 8-9-10     | 3-4.5-6                | 1800        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1055G2 | 55          | 10-11-12   | 3-4.5-6                | 1800        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1255G2 | 55          | 12-13.5-15 | 3-4.5-6                | 1800        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1555G2 | 55          | 15-16.5-18 | 3-4.5-6                | 1800        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1855G2 | 55          | 18-19-20   | 3-4.5-6                | 1800        | 2.4        | Wire-guided | 2.8                          |          |
| KZ2055G2 | 55          | 20-21-22   | 3-4-5                  | 1800        | 2.4        | Wire-guided | 3.2                          |          |
| KZ0655H2 | 55          | 6-7-8      | 3-4.5-6                | 2400        | 2.4        | Wire-guided | 2.8                          |          |
| KZ0855H2 | 55          | 8-9-10     | 3-4.5-6                | 2400        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1055H2 | 55          | 10-11-12   | 3-4.5-6                | 2400        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1255H2 | 55          | 12-13.5-15 | 3-4.5-6                | 2400        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1555H2 | 55          | 15-16.5-18 | 3-4.5-6                | 2400        | 2.4        | Wire-guided | 2.8                          |          |
| KZ1855H2 | 55          | 18-19-20   | 3-4.5-6                | 2400        | 2.4        | Wire-guided | 2.8                          |          |
| KZ2055H2 | 55          | 20-21-22   | 3-4-5                  | 2400        | 2.4        | Wire-guided | 3.2                          |          |